



PHOENIX

MICROWAVE MUFFLE FURNACE
Fast, Accurate Ash Analysis



- a** Heating element placed in furnace walls for optimum heating efficiency
- b** Dual protection system prevents microwave door from closing without the furnace door in place
- c** Large display for easy viewing
- d** Programmable temperature control allows simple single set point temperature or stored multiple stage programs
- e** NIST-traceable dual thermocouple allows rapid verification and calibration of the furnace temperature for ISO/GLP procedures
- f** High speed exhaust eliminates need for unit placement in a fumehood
- g** Furnace temperature uniformity exceeds requirements of standard methods (USP, ASTM, AOAC)
- h** Use any kind of crucible, including metal!
- i** Store up to 20 methods
- j** IR sensor shuts down microwaves in the event of overheating
- k** Provides safe storage for furnace door
Protects operators from hot inner surface

Two furnace configurations available!



High Temperature
1200 °C Furnace holds up to 8 (25-mL) crucibles



High Capacity
1000 °C Furnace Holds up to 15 (25-mL) crucibles Ideal for higher-throughput laboratories

PHOENIX FEATURES

Fast & Efficient

Phoenix Muffle Furnaces perform many high temperature applications up to 97% faster than traditional muffle furnaces, giving you more time to make adjustments to your process and reduce out-of-specification product.

Accurate

Phoenix furnaces satisfy standard methods that require electrically heated furnaces. They have built-in calibration software and NIST traceable accessories are available.

Safe

Phoenix reduces exposure to fumes and heat with a built-in exhaust system that does not require placement of the unit in a fume hood. These rugged, durable furnaces feature door interlocks and built-in system diagnostics.

Easy-to-Use

Auto-start software allows pre-programmable warm-up and shut down. Phoenix furnaces store up to 20 user-programmed methods.

- Up to 97% faster than conventional muffle furnaces
- Accurate results in minutes
- 2 furnace sizes to choose from
- Clean, cool operation
- Rugged and durable
- Reduces exposure to fumes and heat
- Quartz fiber crucibles allow rapid cooling and eliminates risks of burns from handling
- Platinum, porcelain, and graphite crucibles allowed
- Airwave version available
- Sulfated Ashing option available

QUICK & EASY STEPS

ACCURATE ASH ANALYSIS Doesn't Get Any Easier Than This

1



Select the ashing program for your sample using the keypad on the Phoenix Ashing Furnace

2



Weigh your sample into a crucible

3



Place samples in furnace

4



Press "Start"



Over 50,000 systems sold worldwide



CEM has been an ISO-certified facility since 1994



All systems serviced & supported by experts with an average of 15 years of experience



CEM invests 11% of annual revenue into R&D, the result... 11 R&D 100 awards



IQ/OQ/PQ Validation by certified CEM Technicians

“

Our commitment to you doesn't end when your system is shipped; it begins.

~Michael J. Collins
President & CEO, CEM

”

CEM Corporation:
PO Box 200
Matthews, NC 28106
United States
800-726-3331
704-821-7015
Fax: 704-821-7894
info@cem.com
www.cem.com

France: CEM µWave S.A.S.
Immeuble Ariane
Domaine Technologique de Saclay
4, rue Rene' Razel, 91892 ORSAY Cedex
33 (01) 69 35 57 80 • Fax: 33 (01) 60 19 64 91
info.fr@cem.com

Germany, Austria, Switzerland: CEM GmbH
Carl-Friedrich-Gauss-Str.9, 47475 Kamp-Lintfort
(49) 2842-9644-0 • Fax: (49) 2842-9644-11
info@cem.de

Ireland: CEM Technology (Ireland) Ltd.
Sky Business Centre, 9a Plato Business Park,
Damastown, Dublin 15
+353 (0) 1 885 1752 • Fax: +353 (0) 1 885 1601
info.ireland@cem.com

Italy: CEM S.R.L.
Via Dell' Artigianato, 6/8
24055 Cologno al Serio (Bg)
(39) 35-896224 • Fax: (39) 35-891661
info.srl@cem.com

Japan: CEM Japan K.K.
2-18-10 Takanawa, Minato-ku, Tokyo 108-0074
+81-3-5793-8542 • Fax: +81-3-5793-8543
info@cemjapan.co.jp

United Kingdom: CEM Microwave Technology Ltd.
2 Middle Slade, Buckingham Industrial Estate,
Buckingham MK181WA
(44) 1280-822873 • Fax: (44) 1280-822873
info.uk@cem.com

Worldwide patents issued and pending
©2016 CEM Corporation

Phoenix™ and Phoenix AirWave™ are trademarks of CEM Corporation.

