


Appearance



- Needed setting space is depth 60 cm and width 80 cm.
- PC with Windows XP and USB port is optionally needed

Measurement Example

Measurement and Analyzing software 「PEC-Pro」

■ Reference Measurement Mode

This mode measures data of a reference Si cell. From the measured data and the IPCE spectrum stored in your computer, the software calculates the irradiance at each wavelength.

■ Sample Measurement Mode


This mode measures data of samples, and corrected spectra are showed in real-time on measuring.

■ Automatic calculation of photocurrent density


The photocurrent density of solar cells under the standard solar spectrum (AM1.5G – 100 mW cm⁻²) is automatically calculated from IPCE spectra and showed.

■ Save and load data

Measurement data are easily converted to files in text format, which is loaded by Microsoft Excel and etc. In addition, PEC-Pro software can load text files and process the data.



Each photocurrent density measured by an I-V curve (below) and an IPCE spectrum (top) is compared for discussion.



Example of measurement of a dye-sensitized solar cell using N719 dye.
(top) IPCE spectra measured by PEC-Pro
(bottom) I-V curve measured by PEC-IV2 (sold separately)




Introduction of PEC-S20 IPCE Measurement System

Pecell Technologies, Inc.
Yokohama, Japan

Concept of PEC-S20 IPCE Measurement System

- Measuring IPCE of solar cell **without a shading mask**
A solar cell area is determined not by a shading mask but **by number of photons.**
- Lamp house, power supply, monochromator are compacted in a **unit body**, which is easily placed on a desktop.
- **Contact Irradiation** through an optical fiber
It is easy to calculate a current density without considering actual cell area.
Option: An achromatic lens makes 2mm or 1mm irradiation spot, which is suitable for Small solar cells such as Organic Solar Cells.



Purposes of Measurement of IPCE

- Measure spectral response of photoelectric conversion device
- Estimate a Current density under a standard solar spectrum
In the case of a small cell, difficult to determine an actual cell area, which should be needed in current density calculations.

IPCE spectrum is a plot of a ratio of number of output electrons (current) and input photons (irradiance) against wavelengths.





Figure Examples of IPCE spectrum and solar simulator spectrum

Error of current density calculation without a shading mask (1)




It is very important to consider an actual cell area in current density calculation.

Without a shading mask, calculation of current density is affected by photons from outside of the cell.

Often overestimate the current density, J_{sc} .


Error of Current density without a shading mask (2)

TiO₂ area 5mm x 5mm Dye-sensitized solar cell



| | $J_{sc} / \text{mA cm}^{-2}$ | V_{oc} / V | FF | $Effi. / \%$ |
|---------------------|------------------------------|---------------------|-------------|--------------|
| With mask | 14.0 | 0.79 | 0.74 | 8.20 |
| Without mask | 18.2 | 0.80 | 0.74 | 10.67 |


Contact irradiation through the optical fiber of PEC-S20




All the photons from an optical fiber go into a sample cell.
IPCE is simply calculated by number of photons (irradiance) and electrons (current) against wavelength.
No need to correct actual cell areas.

Measurement of IPCE by PEC-S20

Other system



PEC-S20





PEC-S20

No need to measure actual cell area of solar cells.

Comparison of J_{sc} calculated from IPCE spectrum with I-V curve

$$J_{sc} / mAcm^{-2} = \int \frac{P_{in\lambda} \times \lambda}{1.99 \times 10^{-16}} \times \frac{IPCE_{\lambda}}{100} \times \frac{1}{6.24 \times 10^{18}} d\lambda$$

$J_{sc} = 15.2 \text{ mA cm}^{-2}$ (N719)




J_{sc} can be calculated with using IPCE spectra and standard solar spectrum or solar simulator spectrum.

Achromromatic optics for measurement of IPCE of small solar cells

- For measuring small solar cells with 2 mm or 1 mm Especially for organic solar cells.
 - It is difficult to set appropriately a shading mask on small solar cells. It causes to often overestimate J_{sc} . Sometimes, the error reaches 2 or 3 times higher.
- A solar cell area is determined not by a shading mask but by number of photons. This method makes J_{sc} calculation very easily.



Irradiation to a reference silicon solar cell by the Achromatic optics



A spot of 550 nm light
(Green spot)

How to make shading masks by using an old floppy disk cover

- Shading masks are consumable goods. They are easily made dirty by electrolyte or something.
- A shading mask made by an metal plate is difficult to be prepared by user own.
- A shading mask is possible to be made by what users have. Old black floppy disk drive cover and a stationery punch



The radius of an optical fiber of PEC-S20 is 5.5 mm. Usual stationery punch size is 6 mm or 5.5 mm. A black floppy disk cover is used for mask, and to be a light guide as well as a shading mask

PEC-S20 Models

- **A. Standard**
PEC-S20 (Standard DC measurement)
 - **B. AC Measurement**
 - + Chopper Unit
 - + Lock-in-Amplifier
 - **C. AC Measurement with White bias light**
 - + Chopper Unit
 - + Lock-in-Amplifier
 - + White bias bianntenary optical fiber
 - **D. DC measurement Keithley 2400 model**
 - + Keithley 2400 sourcemeter
 - **E. 1 mm ϕ achromatic lens optics**
 - + Achromatic lens optics
 - **F. 1 mm ϕ achromatic lens optics Keithley 2400**
 - + Achromatic lens optics
 - + Keithley2400
-