



### **FEATURES**

- Ergonomic Italian design for maximum operator
- Communication via Ethernet connection (ModBus TCP/IP) for remote control
- Integrated 10" colour touch panel, simple and
- Software with several account levels able to store recipes and reports
- T/RH% probe already integrated for a completely automatic cycle
- DOP-testable HEPA filter H14 integrated on the unit
- Built-in balance with safety support for H<sub>2</sub>O<sub>2</sub> bottles
- Precision peristaltic pump
- Adjustable blower up to 200 m<sup>3</sup>/h
- Reports exported in pdf format via USB port or printed on remote device

### **TYPICAL APPLICATIONS**

- · Clean room for sterile production
- Isolator for sterility test
- Pass box



# Bioreset Plus

## **OPERATIONS WITH GENERATOR** POSITIONED OUTSIDE THE **APPLICATION**

- Set-up Bioreset Plus outside the application and connect it to the feed-through port
- Turn off air conditioning / HVAC system
- · Secure area with warning signs for personal safety
- Start the recipe on Bioreset Plus
- · The instrument controls and stores relevant cycle parameters of biodecontamination cycle
- Aerate the room/application using the air conditioning / HVAC system
- For personnel safety, upon cycle end check H<sub>2</sub>O<sub>2</sub> gas residue before anybody can enter the room.

### **OPTIONAL ACCESSORIES**

- PPM probe to check H<sub>2</sub>O<sub>2</sub> concentration during the cycle
- Catalyzer for exhaust or return air with internal blower on request
- Special pipe dimensions on request
- Safety devices for personnel safety
- Software compliance with FDA 21 CFR Part 11 regulation
- · Turboflow, a versatile system for ultra-fast H<sub>2</sub>O<sub>2</sub> distribution

Up to 500 m <sup>3</sup>
Up to 200 m³/h, adjustable
Integrated 99,995% MPPS
Integrated 0% 100% RH – step 0,1% -40°C +60°C – step 0,1°C
Integrated 0 g 10.000 g – step 1 g
Integrated 1,5 12 g/min – step 0,5
Optional 0 2.000 PPM
PC or PLC
A ID TODUD
ModBus TCP/IP
Optional Compliance to FDA 21 CFR Part 11 regulation
Optional Compliance to FDA 21 CFR
Optional Compliance to FDA 21 CFR Part 11 regulation
Optional Compliance to FDA 21 CFR Part 11 regulation Up to 4 different computers
Optional Compliance to FDA 21 CFR Part 11 regulation Up to 4 different computers Up to 4 different printers
Optional Compliance to FDA 21 CFR Part 11 regulation Up to 4 different computers Up to 4 different printers 60 mm or tri-clamp 2"1/2 Stainless steel AISI 316
Optional Compliance to FDA 21 CFR Part 11 regulation Up to 4 different computers Up to 4 different printers 60 mm or tri-clamp 2"1/2 Stainless steel AISI 316 and ABS
Optional Compliance to FDA 21 CFR Part 11 regulation  Up to 4 different computers  Up to 4 different printers  60 mm or tri-clamp 2"1/2  Stainless steel AISI 316 and ABS  230V / 50Hz or 110V / 60Hz





