

Taq DNA Polymerase
Technical Manual No. 0213

Update date 08302013

I	Introduction	1
II	Application	1
III	Key Features	1
IV	Ship and Storage	1
V	General PCR Protocol Using <i>Taq</i> DNA Polymerase	1
VI	Order Information	2

I INTRODUCTION

Taq DNA Polymerase is a thermostable DNA Polymerase isolated from an *E. coli* strain that carries the *Taq* DNA polymerase gene. *Taq* DNA polymerase is the most common polymerase used for PCR* reactions

II APPLICATIONS

Taq DNA Polymerase can be used in most applications including the following:

- PCR.
- 3' A-tailing of blunt ends.
- Primer extension.
- DNA sequencing.

III KEY FEATURES

Key features of *Taq* DNA Polymerase:

- Terminal transferase activity. *Taq* DNA Polymerase has terminal transferase activity which results in the addition of a single nucleotide (adenosine) at 3' end of the extension product.
- High-purity. No contamination activity has been detected in standard test reactions.

IV SHIPPING AND STORAGE

This product is shipped on blue ice. Store the product at -20°C.

V GENERAL PCR PROTOCOL USING *Taq* DNA POLYMERASE

1. Thaw all the reagents for PCR on ice. **Vortex to mix to remove concentration gradient** and then spin down briefly.

2. Set up 50 μ l PCR reaction in a thin-wall PCR tube on ice by the following recipe:
 - 5 μ l 10X *Taq* buffer solution containing Mg^{2+} .
 - 1 μ l 10 mM dNTP stock
 - 1 μ l Forward primer (50 μ M)
 - 1 μ l Reverse primer (50 μ M)
 - 2 μ l Template (up to 100 ng/ μ l) sterile or filtered H₂O
 - 39.5 μ l sterile or filtered H₂O
 - 0.5 μ l *Taq* polymerase (5 units/ μ l)
3. Program PCR cyclor as following and start:
 - Initial denaturing: 94°C for 3 minutes
 - Then 30 cycles of:
 - 94°C for 30 seconds
 - 55°C for 45 seconds
 - 72°C for 60 seconds (about 1 kb/minute)
 - Extension: 72°C for 7 minutes
4. When the temperature of PCR cyclor reaches 94°C, put PCR reaction tube in and continue the program.
5. Analyze PCR fragments on a agarose or polyacrylamide gel.

Note:

1. This is a basic protocol. One needs to optimize the reagent concentrations, conditions and parameters.
2. This protocol is for PCR cyclor with a hot lid. Otherwise, mineral oil needs to be added to prevent evaporation.
3. 5% DMSO, 1M betaine, or both can be included in PCR reaction to improve the results when a GC-rich template is used.

VI ORDER INFORMATION

Taq DNA Polymerase, Cat. No. E00007-1000 Cat. No. E00007-50000
Green*Taq* DNA Polymerase, Cat. No. E00043

* The PCR process is covered by U. S. Patent numbers 4683195 and 4683202 issued to Cetus and owned by Hoffman-La Roche Inc. GenScript does not encourage or support the unauthorized use of the PCR process. Use of this product is recommended for persons that either have a license to perform PCR or are not required to obtain a license.

GenScript USA Inc.

860 Centennial Ave., Piscataway, NJ 08854

Tel: 732-885-9188, 732-885-9688

Fax: 732-210-0262, 732-885-5878

E-mail: info@genscript.com

Web: <http://www.genscript.com>

For Research Use Only.

Taq DNA Polymerase without Mg²⁺

Cat. No.: E00008

Size: 1,000 U

Description:

Taq DNA Polymerase is a thermostable DNA Polymerase isolated from an *E. coli* strain that carries the *Taq* DNA polymerase gene. *Taq* DNA Polymerase is the most common polymerase used for PCR* reactions.

Key Feature:

- Terminal transferase activity. *Taq* DNA Polymerase has terminal transferase activity which results in the addition of a single nucleotide (adenosine) at 3' end of the extension product.
- High-purity. No contamination activity has been detected in standard test reactions.

Unit Definition:


One unit is the amount of enzyme that can incorporate 10 nmol of dNTP into acid-insoluble material in 30 minutes at 74°C.

10 X reaction Buffer (without Mg²⁺):

500 mM KCl, 100 mM Tris HCl (pH 9.0 at 25°C), 1% Triton X-100 Buffer. This buffer is optimized for use with 200 μM dNTPs.

Example

PCR performance, activity, nuclease.



Lane	Taq	United Used
1	Leading Brand A	0.1
2		0.25
3		0.5
4	Leading Brand B	0.1
5		0.25
6		0.5
7	GenScript	0.1
8		0.25
9		0.5

Note: If the reaction is performed without this buffer, then add 0.1% Triton X-100 (final concentration) to ensure high activity.

Storage Buffer and Concentration:

The enzyme is delivered in 5 units/μl in 20 mM Tris HCl (pH 8.0), 0.1 mM EDTA, 1 mM DTT, 0.1% Triton X-100 and 50% glycerol.

Storage:

This product can be stored at -20°C for future using.

Formulation:

GenScript *Taq* DNA Polymerase has been formulated using GenScript's proprietary technology. The enzyme can be shipped at room temperature or even 37°C for seven days without any loss of activity.

Applications:

Taq DNA Polymerase can be used in most applications including the following:

- PCR*
- 3' A-tailing of blunt ends
- Primer extension
- DNA sequencin

* The PCR process is covered by US. Patents Nos. 4683195 and 4683202, issued to Cetus and owned by Hoffman-La Roche Inc. GenScript does not encourage or support the unauthorized use of the PCR process. Use of this product is recommended for persons who either have a license to perform PCR or are not required to obtain a license. Sale of this product is restricted to regions or countries where native *Taq* DNA polymerase patents have been invalidated.

Taq DNA Polymerase, concentrated

Cat. No. E00012

Size: 2500 U

Description:

Taq DNA Polymerase is a thermostable DNA polymerase isolated from an *E. coli* strain that carries the *Taq* DNA polymerase gene. *Taq* DNA Polymerase is the most common polymerase used in PCR*. In some cases, such as RAPD PCR, adding large volume of general *Taq* DNA polymerase (5 U/μl), which has a high concentration of glycerol in its storage buffer, will increase the glycerol concentration in the reaction mix, interfering with PCR performance. The use of concentrated *Taq* DNA Polymerase (25 U/μl), with a far slimmer dose of glycerol, can prevent poor PCR efficiency.

Note: Concentrated *Taq* DNA Polymerase (GenScript, E00012) is supplied with 10X reaction buffer containing 15 mM magnesium chloride. The dNTP (10 mM) mixture may be ordered separately (See related products).

Key Feature:

- Terminal Transferase Activity: A single nucleotide (adenosine) is added to the 3' end of the extension product.
- High-Purity: No contamination activity has been detected in standard test reactions.

Concentration:

Supplied in 25 units/μl in 20 mM Tris HCl (pH 8.0), 0.1 mM EDTA, 1 mM DTT, **0.1% Triton X-100** and 50% glycerol.

Unit Definition:

One unit is defined as the amount of enzyme that can incorporate 10 nmol of dNTP into acid-insoluble material in 30 minutes at 74°C.

10 X reaction Buffer (with Mg²⁺)

500 mM KCl, 100 mM Tris HCl (pH 9.0 at 25B0C), 15 mM MgCl₂, 1% Triton X-100 Buffer. This buffer is optimized for use with 200 μM dNTPs.

Important:

If another reaction buffers are used with *Taq* DNA Polymerase, Triton X-100 must be added to a final concentration of 0.1% to ensure high enzyme activity with *Taq* DNA Polymerase;

Storage:

Store the product at -20°C.

Formulation:

GenScript *Taq* DNA Polymerase has been formulated using GenScript's proprietary technology. The enzyme can be shipped at room temperature or stored at 37°C for seven days without any significant loss of activity.


Applications:

The applications of *Taq* DNA Polymerase are as follows:

- PCR*
- 3' A-tailing of blunt ends (T/A-cloning)
- Primer extension
- DNA labeling reactions

QC Tests

PCR performance, activity, nuclease.



Lane	Taq	Unit Used
1		0.1
2	Leading Brand A	0.25
3		0.5
4		0.1
5	Leading Brand B	0.25
6		0.5
7		0.1
8	GenScript	0.25
9		0.5

* The PCR process is covered by US. Patents Nos. 4683195 and 4683202, issued to Cetus and owned by Hoffman-La Roche Inc. GenScript does not encourage or support the unauthorized use of the PCR process. Use of this product is recommended for persons who either have a license to perform PCR or are not required to obtain a license. Sale of this product is restricted to regions or countries where native *Taq* DNA polymerase patents have been invalidated.