

**Technical Note** 

## Nicomp DLS Intensity-, Volume-, and Number-weighted Gaussian Distributions

The Nicomp Model 380 measures particle size by detecting the scattered light (at a fixed angle) produced by an ensemble of particles suspended in solvent. The Intensity  $I_s(t)$  fluctuates in time due to diffusion of the particles. There is a well defined characteristic lifetime of the *intensity fluctuations*, which is inversely proportional to the particle diffusivity. We computer the autocorrelation function C(t') of the fluctuating intensity  $I_s(t)$ , obtaining a decaying exponential curve in time. From the decay time constant  $\tau$ , we obtain the particle diffusivity, D. Suing the Stokes-Einstein relation:

## D=kT/6πηR

where k = Boltzmann's constant  $\eta$  = viscosity of solvent (viscosity of water is 0.933 x 10<sup>-2</sup> poise at 23°C) T = temperature (°K = °C + 273)

We then compute the radius, R, assuming a sphere.

Since the autocorrelation function C(t') is constructed from the original scattered intensity value as a function of time, the quadratic fit and the corresponding Gaussian-like representation of the distribution of particle diffusivities (and, ultimately, diameters) reflect the fact that the D (diffusivities) contributions, or R (radius) contributions, are weighted by their corresponding scattering intensities.

There are six pieces of quantitative information contained in the summary display (Figure 1) of the Gaussian Analysis:

Intensity-weighted Mean: nm	Mean Diameter: 219
Width of the Distribution:	Standard Deviation: 34.2 nm
Goodness of Fit: (Qualitative amount)	Gaussian X <sup>2</sup> Chi Squared: 0.55
Aggregates of Dirt:	Auto Baseline Adjustment: 0.538%

## Actual Diffusion Coefficient: Mean Diffusion

Amount of Data Taken:

Data

Particle S	izing Systems	Title: Liposyn 10% f Recipe: Database: Import_Ddata.	at emulsion adt		
User: Run Date/Time: 2/11/2010, 12:53 Report Date/Time: 2/11/2010, 1:22: Report ID: 133	3:54 P 26 PM	Model number: Serial number: Manufacture date: Solid Particle Ana	alysis		
Comment: File Import from F:\\Fa	at_emulsion.399 Original Dat	eStamp: 11:22:50 2/3/2010	D		
	:	Summary			
Runtime: 0.7:57 Channel 1 data: 769102 (x1000) K Average intensity: 278681 kHz Sensitivity: 0		Fit error: 4.263 µm Residual: 7.370 Gaussian ¥ (Chi squared): 0.55 Auto baseline adjustment: 0.538 %			
Gaussian distribution data	NiComp dist	ribution data			
Baussian distribution data Intensity weighted Mean diameter: 219.0 nm	NiComp dist	ibution data Peak 1 data ity weighted Diameter: 219.1 nm	Peak 2 data 0.0 nm	Peak 3 data 0.0 nm	
Gaussian distribution data Intensity weighted Mean diameter: 219.0 nm S.D. (nm/%): 34.2 nm	NiComp dist Intens (15.6 % %)	ibution data Peak 1 data Ity weighted Diameter: 219.1 nm Percent: 100.00%	Peak 2 data 0.0 nm 0.00%	Peak 3 data 0.0 nm 0.00%	
Saussian distribution data Intensity weighted Mean diameter: 219.0 nm S.D. (nm%): 34.2 nm i Volume weighted	NiComp dist Intens (15.6 % %) Volui	ribution data Peak 1 data Ity weighted Diameter: 219.1 nm Percent: 100.00% ne weighted	Peak 2 data 0.0 nm 0.00%	Peak 3 data 0.0 nm 0.00%	
Gaussian distribution data Intensity weighted Mean diameter, 219,0 nm S.D. (nm%): 34.2 nm : Volume weighted Mean diameter. 213.4 nm S.D. (nm%): 33.3 nm ;	NiComp dist Intens (15.6 % %) Volur (15.6 % %)	ribution data Peak 1 data Diameter: 219.1 nm Percent: 100.00% ne weighted Diameter: 206.8 nm Percent: 100.00%	Peak 2 data 0.0 nm 0.00% 0.0 nm 0.00%	Peak 3 data 0.0 nm 0.00% 0.0 nm 0.00%	
Gaussian distribution data Intensity weighted Mean diameter: 219.0 nm S.D. (nm?%): 34.2 nm Volume weighted Mean diameter: 213.4 nm S.D. (nm?%): 33.3 nm Dumber weighted	NiComp dist Intens (15.6 % %) (15.6 % %)	ribution data Peak 1 data lity weighted Diameter: 219.1 nm Percent: 100.00% ne weighted Diameter: 206.8 nm Percent: 100.00%	Peak 2 data 0.0 nm 0.00% 0.0 nm 0.00%	Peak 3 data 0.0 nm 0.00% 0.0 nm 0.00%	
Gaussian distribution data Intensity weighted Mean diameter: 219.0 nm S.D. (nm%): 34.2 nm Volume weighted Mean diameter: 213.4 nm S.D. (nm%): 33.3 nm Number weighted Mean diameter: 196.7 nm S.D. (nm%): 30.7 nm	NiComp dist Intens (15.6 % %) (15.6 % %) (15.6 % %)	ribution data Peak 1 data Ity weighted Diameter: 219.1 nm Percent: 100.00% Diameter: 206.8 nm Percent: 100.00% Diameter: 184.0 nm Percent: 100.00%	Peak 2 data 0.0 nm 0.00% 0.0 nm 0.00%	Peak 3 data 0.0 nm 0.00% 0.0 nm 0.00%	
Gaussian distribution data Intensity weighted Mean dameter: 719 0 nm S.D. (mm%): 342 nm Volume weighted Mean diameter: 7134 nm Number weighted Mean diameter: 166 7 nm S.D. (mm%): 30.7 nm	NiComp dist (15.6 % %) (15.6 % %) (15.6 % %) Numl	Ibution data Peak 1 data Ity weighted Diameter: 219.1 nm Percent: 100.00% ne weighted Diameter: 206.8 nm Percent: 100.00% ber weighted Diameter: 184.0 nm Percent: 100.00%	Peak 2 data 0.0 nm 0.00% 0.0 nm 0.00%	Peak 3 data 0.0 nm 0.00% 0.00% 0.0 nm 0.00%	

Figure 1: Summary of the Gaussian Distribution



Figure 2: Intensity-weighted Gaussian Distribution



The problem is the value of the Mean Diameter is intensity-weighted. Most users do not care which particle number-weighted distribution: contributed the most scattering, intensity. What they would like to know is what size contributes the most volume, volume-weighted or where the largest population of particles is located, number-weighted.

First, let us look at the simplest case: the Raleigh region. This is defined where the particle diameter is much less than  $\tau$ , the wavelength of light (typically < 150 nm). The particles scatter light as point scatterers (equal in all directions).

Let us determine the behavior between intensity-, volume -, and number-weighted by reviewing the relationship between particle size and scattering intensity. For particle sizes much smaller that the wavelength of light (< 150 nm) you are in the Raleigh region where:

 $f(n_p n_s) x (Mw)^2 x I_0 or$  $I_s = g(n_p n_s) \times V^2 \times I_o$  $n_p$  = is the index of refraction of the particle where  $n_s$  = is the index of refraction of the solvent  $I_0$  = laser intensity

And where the contribution C(t') of particle decaying exponential, corresponding to a given diameter, d<sub>i</sub>, should be weighted by the factor  $N_i(Vi)^2$  where  $N_i$  is the number of particles having diameter,  $d^i$ , and  $V_i$  is their individual volume, which is proportional to (di)<sup>3</sup>.

Hence, each diameter "slice" of the intensity-weighted distribution (Figure 2) has associated with it this factor  $N_i d_i^6$ .

Intensity-weighted  $\alpha N_i V_i^2$  or  $N_i d_i^6$ 

To convert this to a volume-weighted (Figure 3) distribution we remove a factor of  $V_i$  or  $d^3$  from each slice of the distribution.

Volume-weighted  $\alpha N_i V_i^3$ 



Figure 3: Volume-weighted Gaussian Distribution

We perform one additional division of Vi to obtain the

Number-weighted a Ni

For sizes larger than 150 nm, the rules of Mie scattering take effect when intra-particle interference effects become important.



Figure 4: Number-weighted Gaussian Distribution

The Nicomp uses the "rules" of light scattering making corrections for Mie scattering. In order to obtain the volume- and number-weighted plot distributions from the intensity-weighted plot the relationships are as follows:

> Intensity-weighted  $\alpha N_i V_i^2 G_i$ Volume-weighted a NiViGi Number-weighted α N<sub>i</sub>G<sub>i</sub>

where G is a factor that varies between unity and zero depending on the particle size, laser wavelength, and angle of detection.

That is why there are some cases in larger particle size distributions where the volume-weighted mean diameter is larger than the intensity weighted mean.

Also keep in mind, the amount that the mean diameter of the distribution shifts depends on the width of the intensity-weighted distribution. An extreme case is where you have a mono-sized peak; the intensity-, volume- and number-weighted means would all be the same.

## **Particle Sizing Systems**

8203 Kristel Circle, New Port Richey, FL 34668 Phone: 727•846•0866 | Fax: 727•846•0865 Website: www.pssnicomp.com E-mail: sales@pssnicomp.com

Particle Sizing Systems



Building solutions one particle at a time