

A

Sterl IT

< > €





# Sterilizing by IR-HotSpot



111



New







Efficient use of energy



No air turbulences

## Infrared loop sterilizer - using IR-HotSpot technology

#### High-Speed annealing and sterilizing the smartest way

Are you looking for a fast, efficient and flameless product to anneal or sterilize your inoculation loops and delicate micro instruments? Then our electrical *SteriMax basic / smart* should be your first choice – ideally suited for use in any type of lab, anaerobic work environment and Laminar Flow Cabinet!

#### Smart infrared light - instantly ready to work

Specifically focused infrared light generates an IR-HotSpot in which your inoculation loop is sterilized at temperatures from 750 to 1000°C after only 5 to 10 seconds. The *SteriMax basic / smart* just needs to be plugged in and will be immediately ready to use, without any warm-up period!

Our highly sophisticated annealing tube is based on many years of experience. It is made of special quartz glass. The tube is closed at one end preventing any risk of contamination with pathogenes. The annealing tube can be slid out in a matter of seconds for easy cleaning.

#### Smart operation - comfortable working

*SteriMax basic / smart* allows you to fully focus on your core task. Hardly any attention is needed for inserting your inoculation loop into the annealing tube – which is easily accessible from the front of the flat glass panel. Depending on the model you can start the sterilization process via touch panel, IR-Sensor or foot pedal\*.

#### Smart use of energy - minimum heat dissipation

"Thermocontrol" is the key element for an efficient use of energy. It allows a quick but smooth softstart of the infrared light, a short heating period and the re-use of residual heat for the next sterilization process. Whenever you pause, your energy consumption and heat dissipation will pause as well.

### Smart design – outstanding features

*SteriMax smart / basic* is presented in an elegant stainless steel housing with a safety glass touch panel designed for extreme laboratory environments with many unique features.



- Instantly ready to work without any warm-up period
- No ranger of combustion on the housing (Cool-Touch)
- Comfortable touch operation
- Efficient use of energy through Thermocontrol
- Annealing tube made of special quartz glass
- Very simple cleaning of the device and annealing tube





The sterilization process starts with just a fingertip on the touch panel or by pushing the foot pedal\*.

A fixed timer (10 sec.) controls the annealing of your instruments automatically.

Thermocontrol uses the residual heat and automatically shortens the energy supply for the next sterilization.

The course of the sterilization can be followed from start to end by looking upon the inoculation loop symbol.

After that your inoculation loop is ready for use!

\*optional

	⊖ <b>Steri</b> Max	<del>smart s <u>Steri</u>Max</del>
Programs:		
Start via IR-Sensor	4	$\checkmark$
Start via touch	$\checkmark$	
Start via foot pedal*	✓	
2 user accounts		$\checkmark$
Sterilization at start-up	$\checkmark$	$\checkmark$
Adjustable sterilization- / cool down period		$\checkmark$
Adjustable sensor range		$\checkmark$
Timer display		~
Sterilization display	$\checkmark$	~
No warm-up period	$\checkmark$	~
Thermocontrol	$\checkmark$	$\checkmark$
Residual heat display	$\checkmark$	$\checkmark$
Buzzer		$\checkmark$
Holding device for 3 WLD-TEC inoculation loop holder		$\checkmark$
Cool-Touch casing	$\checkmark$	$\checkmark$

\*optional



Touchless IR-Sensor technology starts the sterlization process automatically.

The timer can be viewed and adjusted (exact to the second) for 2 different user accounts simply by touching the panel.

Thermocontrol is using the residual heat and automatically reduces the time of sterilization up to 5 seconds.

2 adjustable timers indicate the completion of the sterilization period (5 - 10 sec.) and cooling down process (0 - 25 sec.) via the display and an audible signal.

After that your inoculation loop is ready for use!





#### The range:

Infrared light technology and quartz glass annealing tube Replacement quartz glass annealing tube (smart) Inoculation loop holder with inoculation loop Ø 3 mm Holding device for 3 WLD-TEC inoculation loop holders (smart) Instruction manual, 2-year warranty (except parts of wear and tear)

for 230 V mains voltage for 120 V mains voltage for 100 V mains voltage

Art.-No. 4.101.000 Art.-No. 4.102.000 Art.-No. 4.103.000



for 230 V mains voltage for 120 V mains voltage for 100 V mains voltage Art.-No. 4.001.000 Art.-No. 4.002.000 Art.-No. 4.003.000



made in Germany

CE



## **Technical Data**

Sterimax smart IR-Sensor:

Touch operation:

Automatic unit switch off:

Sterimax basic Touch operation: Foot pedal\* operation: Automatic unit switch off:

All models

Status display: Thermocontrol: Sterilization at start-up:

Annealing tube Quartz glass:

Safety features Residual heat display: Temperature monitor:

Electrical Power supply:

Power consumption: Protection class: Heat source:

Mechanical Quartz glass annealing tube: Casing:

Measurements (w x h x d): Weight: Licences CE:

EU guidelines: Protected Design start of sterilization sensor range (8-stages adjustable) ON / OFF, 2 user accounts sterilization timer 5 - 10 sec. inoculation loop cooling down timer 0 - 25 sec. with acoustic signal (disengageable) 15, 30 or 60 min.

start of sterilization (10 sec.) start of sterilization (10 sec.) 15 min.

optical display of the sterilization progress dynamic sterilization timer sterilization of the annealing tube before usage

750 °C - 1000 °C sterilization temperature

indicates a hot annealing tube thermal circuit breaker

230 V mains voltage: 220 - 240 V 50/60 Hz 2.6 A 120 V mains voltage: 112 - 128 V 50/60 Hz 5.1 A 100 V mains voltage: 95 - 105 V 50/60 Hz 6 A approx. 600 W / short time (5 - 10 sec.) IP 20 IR-Halogen lamp

outer-Ø 19 mm, inner-Ø 16 mm, length 112 mm heat resistant glass / stainless steel, UV and solvent resistant 110 x 170 x 180 mm approx. 1800 g

EN 61326-1, EN 61010-1, EN 61010-2-010, EN 62471 2014/30/EU, 2014/35/EU, 2011/65/EU, 2006/25/EG D DM/082 011 (smart) U DE 20 2016 101 067 (basic)

\* foot pedal optional



Tray, 3 places stainless steel, for right or left side Art.-No. 4.000.341



Quartz glass annealing tube Art.-No. 4.000.400



Foot pedals only SteriMax basic

(1) Foot pedal stainless steel, cable length 3.0 m Art.-No. 6.000.402

(2) Foot pedal radio set RF-Stick only available in EU countries

Art.-No. 6.000.404-RF (3) Foot pedal Mini

plastic, cable length 3.0 m Art.-No. 6.000.403



Sales Department: Beethovenstr. 3 D-37085 Göttingen Tel: +49 (0)551 / 793789 Fax: +49 (0)551 / 793707

Email: sales@wld-tec.com • Internet: http://www.wld-tec.com (03/18) Subject to technical changes.



Magnetic instrument support stainless steel Art.-No. 4.000.350



Inoculation loop holder stainless steel, with sleeve nut, for inoculation loops  $\emptyset$  0.6 - 1 mm

length 245 mm Art.-No. 6.000.360 length 215 mm Art.-No. 6.000.365



Inoculation loops, wire Ø 0.6 mm

Special stainless steel Ø 3 mm Art.-No. 6.000.373 Ø 5 mm Art.-No. 6.000.375

Platinum / Iridium Ø 1 mm Art.-No. 6.000.381 Ø 3 mm Art.-No. 6.000.383 Ø 5 mm Art.-No. 6.000.385

 Production & Service:

 Halle-Kasseler-Straße 49

 D-37318 Arenshausen

 Tel:
 +49 (0)36081 / 68940

 Fax:
 +49 (0)36081 / 68942