

Phone: (717) 261-9700
Fax: (717) 261-1730

1671 Orchard Drive
Chambersburg, PA 17201

[Home](#)[About Us](#)[The TBJ Way](#)[Products](#)[Contact Us](#)**TBJ**
INCORPORATED[< Back](#)

30-60-S DD/BD Dual Draft - Fast Facts:

- The 30-60-S DD/BD Dual Draft Station is designed for gross dissection, trimming and necropsy procedures. It employs both downdraft and backdraft ventilation. It is shown here with dual task lights, safety splash shield, formalin cup sink with collection carboy, magnetic tool bar and formalin dispensing carboy.
- Overall dimensions: 30" wide x 60" x 34-40" high work surface
- 41" wide x 24" downdraft work area
- 41" wide x 10" high adjustable backdraft grille
- Spray gun, flush manifold, sink with hot/cold water faucet, fluorescent light fixture standard
- 4" x 40" duct connection



Model 30-60-S DD/BD Dual Draft Workstation



Technical Specifications



Copyright © 2008 TBJ Incorporated

Site Design by Global Data Consultants

[Terms & Conditions](#)

