



## multiSUB<sup>®</sup> Choice

### FEATURES:

- Three tray options
- Run up to 210 samples
- Low buffer volumes
- Multichannel pipette compatible combs for speed loading

The multiSUB<sup>®</sup> Choice offers a wide degree of versatility. Three tray options are available – 15 x 7cm, 15 x 10cm and 15 x 15cm – allowing the choice of one, two or all three gel length options at the time of purchase. Further purchases of additional accessories are no longer required. Maximising comb and tray options allow up to 210 samples to be resolved per gel. The 15cm total run length allows restriction fragment or other close MW sample bands to be easily

separated and identified. Speed loading is accomplished using 10, 14, 16, 18, 28 or 30 sample multi-channel pipette compatible combs. multiSUB<sup>®</sup> choice stretch units are available with optional 15 x 20cm and 15 x 25cm gel trays and four 28-sample combs for those researchers wanting to perform higher resolution separation of more samples over a longer distance. multiSUB<sup>®</sup> Choice Trio includes all 3 tray sizes for optimum versatility and value.

## MULTISUB® CHOICE COMBS

Code	Description	Sample volume for a 5mm thick gel
MS15-1-0.75	Comb Prep 1, Marker 1, 0.75mm thick	371µl
MS15-2-0.75	Comb Prep 2, Marker 2, 0.75mm thick	169µl
MS15-4-0.75	Comb Prep 4, Marker 2, 0.75mm thick	91µl
MS15-10-0.75	Comb 10 sample, 0.75mm thick	34µl
MS15-12-0.75	Comb 12 sample, 0.75mm thick	30µl
MS15-20-0.75	Comb 20 sample, 0.75mm thick	16µl
MS15-35-0.75	Comb 35 sample, 0.75mm thick	7µl
MS15-10MC-0.75	Comb 10 sample MC, 0.75mm thick	22µl
MS15-14MC-0.75	Comb 14 sample MC, 0.75mm thick	22µl
MS15-16MC-0.75	Comb 16 sample MC, 0.75mm thick	20µl
MS15-18MC-0.75	Comb 18 sample MC, 0.75mm thick	8µl
MS15-28MC-0.75	Comb 28 sample MC, 0.75mm thick	8µl
MS15-30MC-0.75	Comb 30 sample MC, 0.75mm thick	9µl
MS15-1-1	Comb Prep 1, Marker 1, 1mm thick	495µl
MS15-2-1	Comb Prep 2, Marker 2, 1mm thick	225µl
MS15-4-1	Comb Prep 4, Marker 2, 1mm thick	122µl
MS15-10-1	Comb 10 sample, 1mm thick	45µl
MS15-12-1	Comb 12 sample, 1mm thick	41µl
MS15-20-1	Comb 20 sample, 1mm thick	21µl
MS15-35-1	Comb 35 sample, 1mm thick	10µl
MS15-10MC-1	Comb 10 sample MC, 1mm thick	29µl
MS15-14MC-1	Comb 14 sample MC, 1mm thick	29µl
MS15-16MC-1	Comb 16 sample MC, 1mm thick	27µl
MS15-18MC-1	Comb 18 sample MC, 1mm thick	11µl
MS15-28MC-1	Comb 28 sample MC, 1mm thick	11µl
MS15-30MC-1	Comb 30 sample MC, 1mm thick	13µl
MS15-1-1.5	Comb Prep 1, Marker 1, 1.5mm thick	743µl
MS15-2-1.5	Comb Prep 2, Marker 2, 1.5mm thick	338µl
MS15-4-1.5	Comb Prep 4, Marker 2, 1.5mm thick	182µl
MS15-10-1.5	Comb 10 sample, 1.5mm thick	68µl
MS15-12-1.5	Comb 12 sample, 1.5mm thick	61µl
MS15-20-1.5	Comb 20 sample, 1.5mm thick	32µl
MS15-35-1.5	Comb 35 sample, 1.5mm thick	15µl
MS15-10MC-1.5	Comb 10 sample MC, 1.5mm thick	44µl
MS15-14MC-1.5	Comb 14 sample MC, 1.5mm thick	44µl
MS15-16MC-1.5	Comb 16 sample MC, 1.5mm thick	41µl
MS15-18MC-1.5	Comb 18 sample MC, 1.5mm thick	17µl
MS15-28MC-1.5	Comb 28 sample MC, 1.5mm thick	17µl
MS15-30MC-1.5	Comb 30 sample MC, 1.5mm thick	19µl
MS15-1-2	Comb Prep 1, Marker 1, 2mm thick	990µl
MS15-2-2	Comb Prep 2, Marker 2, 2mm thick	450µl
MS15-4-2	Comb Prep 4, Marker 2, 2mm thick	243µl
MS15-10-2	Comb 10 sample, 2mm thick	90µl
MS15-12-2	Comb 12 sample, 2mm thick	81µl
MS15-20-2	Comb 20 sample, 2mm thick	43µl
MS15-35-2	Comb 35 sample, 2mm thick	20µl
MS15-10MC-2	Comb 10 sample MC, 2mm thick	59µl
MS15-14MC-2	Comb 14 sample MC, 2mm thick	59µl
MS15-16MC-2	Comb 16 sample MC, 2mm thick	54µl
MS15-18MC-2	Comb 18 sample MC, 2mm thick	23µl
MS15-28MC-2	Comb 28 sample MC, 2mm thick	23µl
MS15-30MC-2	Comb 30 sample MC, 2mm thick	25µl

## ORDERING INFORMATION


MSCHOICE7 multiSUB® Choice	15 x 7cm UV Tray, 2 x 20 sample combs, loading guides and dams
MSCHOICE10 multiSUB® Choice	15 x 10cm UV Tray, 2 x 20 sample combs, loading guides and dams
MSCHOICE15 multiSUB® Choice	15 x 15cm UV Tray, 2 x 20 sample combs, loading guides and dams
MSCHOICETRIO multiSUB® Choice Trio	15 x 7, 10 & 15cm UV Tray, 2 x 20 sample combs, loading guides and dams
MSCHOICETRIO15 multiSUB® Choice Trio 15	3 x 15cm UV tray, 2x 20 sample combs, loading guides and dams
MSCHOICEST20 multiSUB® Choice Stretch 20	15 x 20cm UV tray, 4 x 28 sample combs, loading guides and dams
MSCHOICEST25 multiSUB® Choice Stretch 25	15 x 25cm UV tray, 4 x 28 sample combs, loading guides and dams

## TECHNICAL SPECIFICATIONS


Gel dimensions (w x l)	15 x 7cm    15 x 10cm 15 x 15cm    15 x 20cm    15 x 25cm
Unit dimensions (w x l x h)	17.5 x 26.5 x 9cm 17.5 x 41 x 9cm (Stretch models only)
Max. sample capacity	15 x 7cm tray – 70 samples 15 x 10cm tray – 140 samples 15 x 15cm tray – 210 samples 15 x 20cm tray – 280 samples 15 x 25cm tray – 350 samples
Buffer volume	500ml (choice) 1000ml (choice stretch)
Combs available: No. of samples	1, 2, 4, 10, 10MC, 12, 14MC, 16, 18MC, 20, 28MC, 30MC, 35
Thicknesses	0.75, 1, 1.5, 2mm

## MULTISUB® CHOICE ACCESSORIES


MS15-UV7	15 x 7cm UV Tray	MS15-UVDAM	Casting Dams
MS15-UV10	15 x 10cm UV Tray	MS15-PE	Positive Electrode
MS15-UV15	15 x 15cm UV Tray	MS15-NE	Negative Electrode
MS15-UVST20	15 x 20 UV Tray	MS15-LG	Adhesive Loading Guides
MS15-UVST25	15 x 25 cm V Tray	MS15-WP	Viewing Platform
MS15-UVS	15cm UV Gel Scoop	MSCHOICECP	Cool-Pack and Platform
MS15/20-FC	multiSUB® Choice/Maxi Flexi caster		



Condensation free lids available



multiSUB® Choice stretch version



## MULTISUB® CHOICE ACCESSORIES



Clever Scientific Ltd  
Unit 4 Triton Park, Brownsover Road  
Swift Valley, Rugby, Warwickshire CV21 1SG  
United Kingdom  
t: +44 (0)1788 565300 f: +44 (0)1788 552822  
info@cleverscientific.com www.cleverscientific.co.uk

Distributor