



# Media Formulation

## 0.25% Trypsin-EDTA, 1X Cat. No.: CC508-0100

Components	Concentration (mg/L)
<b>Inorganic Salts</b>	
Potassium Chloride (KCl)	400.0
Potassium Phosphate monobasic (KH <sub>2</sub> PO <sub>4</sub> )	60.0
Sodium Bicarbonate (NaHCO <sub>3</sub> )	350.0
Sodium Chloride (NaCl)	8000.0
Sodium Phosphate dibasic (Na <sub>2</sub> HPO <sub>4</sub> ·7H <sub>2</sub> O)	90.0
<b>Other Components</b>	
D-Glucose (Dextrose)	1000.0
EDTA·4Na·2H <sub>2</sub> O	380.0
Phenol Red	10.0
Trypsin	2500.0