



Media Formulation

TrypRC™ Clear, 1X

Cat. no.: CC512-0100

Inorganic Salts (g/liter)	Concentration (g/L)
KCl	0.2
Potassium Phosphate monobasic (KH_2PO_4)	0.2
Sodium Chloride (NaCl)	8
Sodium Phosphate dibasic ($\text{Na}_2\text{HPO}_4 \cdot 7\text{H}_2\text{O}$)	2.16

Other (g/liter)	Concentration (g/L)
EDTA	0.4576
rProtease	confidential