

Certificate of Analysis



7.5% Sodium Bicarbonate Solution

Catalog Number: CC530-0100

Lot Number: GL20177

Size: 100 mL

Minimum Shelf Life: Jun-2022

Storage: 2 ~ 30°C

Test	Specifications	Results
Endotoxic Assay	Check & Record	<0.04 EU/ mL
Osmolality	Check & Record	1451
pH	Check & Record	8.2
Sterility Testing	Negative	Negative

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.

Yi-Hsing, Yao
Inspector, Quality
Issued on 2022/12/22