## Specifications: DCT Pads

Rating Voltage:	7.4V	
Battery Capacity:	2200mAh	
Charging Voltage:	8.4 to 8.6V	
Charging Amperage:	1.0 to 1.2A	
Output Power:	15 watt maximum	
Output Power Efficiency:	>98%	

	DCT-15	DCT-20	DCT-25
Dimensions:	15.2 cm W x 20.3 cm L	20.3 cm W x 25.4cm L	35.6 cm W x 35.6 cm L
	(6" x 8")	(8" x 10")	(14" x 14")

| Print | Close Window