



Space Saving

Self-inverted design enables it stand on the bench alone



UV Resistant Housing

Decontamination by UV radiation



Light Weight Only 180g



Easier Access to Filter

Fast removable design of the nose cone



High Quality Battery Indicator

3 colors show the battery status



Reliable Connection

Mini USB-plug charger adapter

Specifications

Pipette Compatibility	1 - 100 ml plastic/glass pipette
Max Aspirating and Dispensing Speed	>7.5 ml/sec
Dispensing Mode	Blow-out/Gravity
Speed Control Wheel	Continuous; +/- indicates Fast/Slow
Number of Dispenses	Approx. 4,500 times (with a 25 ml pipette)
Charging Time	Approx. 4 hours
Charger Adaptor	Input: AC 100 - 240V 50/60 Hz Selectable power plug: US/EU/UK/AU Output: DC 5V/2.0 A
Power Jack	Mini-USB
Filter	Hydrophobic 25 mm diameter, 0.2 μm
Battery	Lithium-ion polymer battery 3.7 V/1,100 mAh
Autoclavable	Nose cone, pipette holder and filter only
Operating Temperature	5 - 40°C
Operating Humidity	RH: 0 - 85%
Weight	180 g, including a battery
Certification	Complies with CE, EN55011, EN61326-2-6 (Class A), EN61010-1 RoHS compliant

Ordering Information

Order No.	Description
BS01-1000	BlueSwan pipette controller, US plug
BS01-2000	BlueSwan pipette controller, EU plug
BS01-3000	BlueSwan pipette controller, UK plug
BS01-4000	BlueSwan pipette controller, AU plug



4F., No. 31, Sec. 2, Chang'an E. Rd.,
Zhongshan Dist., Taipei City 10456, Taiwan (R.O.C.)
Tel: +886-2-8219-1206 Fax: +886-2-8219-1205
Email: sales@blue-raybio.com

