



As a Research-Ready™ solution, Innovive disposable rodent caging offers consistent and standardized housing, pre-bedded cages, and pre-filled bottles – irradiated and delivered to your door.

DISPOSABLE CAGING FOR MICE




INNOCAGE® MOUSE

DESCRIPTION:

Disposable IVC caging for mice, Universal Euro Type II Long
14.7"L x 9.2"W x 5.5"H • 81 in² floor space • sized for 5 mice
100% PET • BPA-FREE • Recyclable

FULL CAGE SETUP:



OPTIONS:

IVC and static cages, empty or pre-bedded bottom, three feeder options, card holder and enrichment

RACKS:

Innorack® (IVC caging only)

WATERING:

Aquavive® water bottle (empty or pre-filled); autowatering and pouches compatible

INNOCAGE® IVC STANDARD LID



DESCRIPTION:

Single filter IVC lid

APPLICATION:

Standard IVC rodent housing; core facilities, quarantine, remote labs, etc.

PRODUCT CODES:

MVX1 (LID+BOTTOM)

MVX1-C8/C8E (LID+PRE-BEDDED BOTTOM, CORN COB OR CORN COB WITH INNORICHMENT™)

MVX1-AD/ADE (LID+PRE-BEDDED BOTTOM, ALPHA-DRI® OR ALPHA-DRI® WITH INNORICHMENT™)



INNOCAGE® IVC LARGE FILTER LID



DESCRIPTION:

Large filter IVC lid with larger ventilation grate and static filter area

APPLICATION:

Facilities that lack redundant power supply

PRODUCT CODES:

MVX6 (LID+BOTTOM)

MVX6-C8/C8E (LID+PRE-BEDDED BOTTOM, CORN COB OR CORN COB WITH INNORICHMENT™)

MVX6-AD/ADE (LID+PRE-BEDDED BOTTOM, ALPHA-DRI® OR ALPHA-DRI® WITH INNORICHMENT™)



INNOCAGE® IVC DUAL FILTER LID



DESCRIPTION:

Dual HEPA filter lid

APPLICATION:

Accommodates research that requires rodent housing with dual filtration

PRODUCT CODES:

MV2-B (LID+BOTTOM)

MV2-B-C8/C8E (LID+PRE-BEDDED BOTTOM, CORN COB OR CORN COB WITH INNORICHMENT™)

MV2-B-AD/ADE (LID+PRE-BEDDED BOTTOM, ALPHA-DRI® OR ALPHA-DRI® WITH INNORICHMENT™)



INNOCAGE® IVC CONTAINMENT LID



DESCRIPTION:

Dual HEPA filter lid with blue cover for easy identification

APPLICATION:

Containment applications

PRODUCT CODES:

MV3 (LID+BOTTOM)

MV3-C8/C8E (LID+PRE-BEDDED BOTTOM, CORN COB OR CORN COB WITH INNORICHMENT™)

MV3-AD/ADE (LID+PRE-BEDDED BOTTOM, ALPHA-DRI® OR ALPHA-DRI® WITH INNORICHMENT™)



IMAGES AND MORE INFORMATION:

WWW.INNOVIVE.COM/MOUSE-CAGE-LIDS.HTML

IMAGES AND MORE INFORMATION:

WWW.INNOVIVE.COM/MOUSE-CAGE-LIDS.HTML

INNOCAGE® STATIC LID



DESCRIPTION:

Static lid with a large filter media and metal grate

APPLICATION:

Static housing, emergency preparedness


PRODUCT CODES:

MS2 (LID+BOTTOM)

MS2-C8/C8E (LID+PRE-BEDDED BOTTOM, CORN COB OR CORN COB WITH INNORICHMENT™)

MS2-AD/ADE (LID+PRE-BEDDED BOTTOM, ALPHA-DRI® OR ALPHA-DRI® WITH INNORICHMENT™)

INNOCAGE® TRANSPORT LID



DESCRIPTION:

Large filter lid without water bottle cavity to accommodate transportation

APPLICATION:

Safe, secure and cost-effective way to transport mice

PRODUCT CODES:

MSX4 (LID+BOTTOM)

MSX4-C8/C8E (LID+PRE-BEDDED BOTTOM, CORN COB OR CORN COB WITH INNORICHMENT™)

MSX4-AD/ADE (LID+PRE-BEDDED BOTTOM, ALPHA-DRI® OR ALPHA-DRI® WITH INNORICHMENT™)

IMAGES AND MORE INFORMATION:

WWW.INNOVIVE.COM/MOUSE-CAGE-LIDS.HTML

CAGE BOTTOM OPTIONS



DESCRIPTION:

100% PET empty universal bottom

PRODUCT CODE: M-BTM



DESCRIPTION:

100% PET universal bottom, pre-bedded with 1/8" corn cob Innorichment™ option available

PRODUCT CODE: M-BTM-C8/C8E



DESCRIPTION:

100% PET universal bottom, pre-bedded with ALPHA-dri® Innorichment™ option available

PRODUCT CODE: M-BTM-AD/ADE

IMAGES AND MORE INFORMATION:

WWW.INNOVIVE.COM/MOUSE-CAGE-BOTTOMS.HTML

MOUSE FEEDERS



DESCRIPTION:

100% PET plastic feeder

PRODUCT CODE: M-FEED-P



DESCRIPTION:

100% PET plastic feeder with metal insert

Stainless steel option available

PRODUCT CODE: M-FEED, M-FEED-SS

IMAGES AND MORE INFORMATION:

WWW.INNOVIVE.COM/MOUSE-FEEDERS.HTML

AQUAVIVE® MOUSE



DESCRIPTION:
300 mL bottle, empty
Amber color option available

PRODUCT CODE: **M-WB-300** (CLEAR),
M-WB-A-300 (AMBER)




DESCRIPTION:
300 mL bottle, pre-filled with
laboratory grade acidified water

PRODUCT CODE: **M-WB-300A**

IMAGES AND MORE INFORMATION:
WWW.INNOVIVE.COM/AQUAVIVE-MOUSE-WATER.HTML

ENRICHMENT & ACCESSORIES



DESCRIPTION:
100% PET reusable cardholder,
accommodates standard cards

PRODUCT CODE: **CRD-HLD-H**



DESCRIPTION:
Innodome™ and Innowheel™,
100% PET, translucent red

PRODUCT CODE: **INNODOME,**
INNOWHEEL



DESCRIPTION:
Innorichment™ nesting sheets,
only available with pre-bedded
corn cob and ALPHA-dri® cages

PRODUCT CODE: **M-BTM-C8E/ADE**

IMAGES AND MORE INFORMATION:
WWW.INNOVIVE.COM/MOUSE-CAGE-ACCESSORIES.HTML



ORDERING YOUR MOUSE CAGING

All Innovive mouse cage
lids are compatible with the
universal Innocage® bottom.

All Innovive mouse caging is
sold in sets. Each set includes
a cage lid and a cage bottom.

***All products delivered
double-bagged and
irradiated***

1. SELECT A LID

- » IVC Single Filter
- » IVC Dual Filter
- » IVC Large Filter
- » IVC Containment
- » Static
- » Transport

2. SELECT BOTTOM OPTION

- » Empty
- » Pre-bedded, 1/8" corn cob
(Innorichment™ optional)
- » Pre-bedded, ALPHA-dri®
(Innorichment™ optional)

3. ADD FEEDER

- » Plastic
- » Plastic w/ metal insert
- » Plastic w/ stainless steel
insert

4. WATER, 300 mL

- » Aquavive®, empty
- » Aquavive®, pre-filled
- » Aquavive® Amber, empty

5. ADD ACCESSORIES

- » Card Holder
- » Innodome™
- » Innowheel™

LID+ BOTTOM SOLD AS A SET