

# SurgiSuite

Multi-Functional Surgical Platform



## Chemical-Resistant Base

*With a stainless steel surgical surface and polycarbonate overlay for added protection*

## Rechargeable Far Infrared Animal Warming Pad

*Recessed in center of platform for a level surgical field*

## Bright LED Lighting\*

*Connects to multi-port USB power strip*

*SurgiSuite shown with optional add-ons  
\*Exact lighting may not be as shown due to product improvements*

Order your SurgiSuite:

**888.572.8887**



**Kent Scientific**  
CORPORATION

## SurgiSuite includes:

- Easy to clean, chemical-resistant surface
- Integrated rechargeable soft warming pad
- Magnetic nose cone/intubation tubing stabilizers, (1-small and 1-large)
- Magnetic limb positioners (4-pack with elastic wires)
- Elastic wire positioner (1-meter)
- Detachable procedure bowl
- Disposable platform covers (5-pack)

1-Year 100% Parts & Labor Warranty

Visit our website for a complete selection of SurgiSuite accessories.

## Specifications:

	SurgiSuite - Standard	SurgiSuite - Large
Base dimensions:	12.90 in x 16.80 in (32.77 cm x 42.67 cm)	15.75 in x 17.80 in (40.01 cm x 45.21 cm)
Stainless steel platform:	11.90 in x 7.90 in (30.23 cm x 20.07 cm)	12.75 in x 10.75 in (32.39 cm x 27.31 cm)
Warming surface:	5.15 in x 4 in (13.08 cm x 10.16 cm)	8.75 in x 6.75 in (22.23 cm x 17.15 cm)

## Optional Add-Ons:



### Monitor and regulate animal temperature

The **PhysioSuite**® with RightTemp® warming and temperature control module easily connects to the SurgiSuite.

[Ask about integrating pulse oximetry.](#)



### Customize your lighting solution\*

LED lighting plugs directly into the included USB power strip.

*\*Exact lighting may not be as shown due to LED lighting updates.*



### Integrate an Intubation Stand into your surgical set-up

Easily intubate mice and rats with an integrated anesthesia mask, neck positioner and fiber optic light.

Order your SurgiSuite today:

**888.572.8887**

[www.kentscientific.com](http://www.kentscientific.com)



**Kent Scientific**  
CORPORATION

04.17.1