

Illuminated Tail Vein Unit

NEW!



PLAS ■ LABS now offers an “Illuminated Rodent Tail Vein Injection Cone”. The new 562-IRC features our “Bright Light” illumination system with rechargeable batteries and recharge unit, clear plastic tail insertion cone, formed white base with suction cup “non-slip” feet, housing batteries, on/off switch, lights, and an electrical plug for recharging the batteries (100-240 VAC capacity.)

What’s Included with the 562-IRC:

- Clear plastic tail insertion cone.
- Formed white plastic base containing batteries, “on / off” switch, and lights.
- Base used for recharging batteries. Two (2) included with system.
- Electrical plug for recharging batteries. (100 to 240 Volt capacity.)

562-IRC

Catalog Number:	562-IRC
Cone Dimensions:	5.75” long x 5” high x 4” opening
Total Footprint Size:	8.5” x 10” (216mm x 254mm)
Light bulb warranty:	30,000 hours
Animal Weight Range:	Mice: 15 to 30 grams Rats: 500 to 600 grams



562-IRC “On”



562-IRC “Off”