

Now available!
Variable speed (VS) model
so you can set aliquoting
speed and volume



ali-Q 2 VS on YouTube

2-in-1 Regular and Repeating Pipet Controller

Faster and Easier Aliquoting. No need to watch the lines. More accurate and repeatable than by eye.

- Push one button to aliquot / multidispense with unmatched repeatability.
- Compatible with all brands and sizes of serological pipets; optimal performance when using Wobble-not™ pipets.
- Innovative system with computer controlled sensors and valves for maximum accuracy and precision.*

*get 2% accuracy and 2% CV at 0.5mL when using Wobble-not™ serological pipets

ali-Q 2 Series Additional Features:

- Now available in Variable Speed (VS) model to accommodate your multidispensing needs.
- In-lab calibration software on-board
- Removable kickstand
- Lighter weight

ali-Q 2 LS (Low speed):

- aliquoting speed ~50% slower than ali-Q 2
- white numbers on dial

ali-Q 2



Removable kickstand



ali-Q 2 VS (Variable speed):

- aliquoting speed adjustable with dial on top

ali-Q 2 series Next generation of ali-Q

Description	Product #	Volume Range
ali-Q 2	2100-0005	0.5 - 5.0 mL
ali-Q 2 LS	2100-0003	0.3 - 3.0 mL
ali-Q 2 VS	2100-1005	0.5 - 5.0 mL

Pipetting is *FASTER* with

The **ONLY** way to perform
Repeat Pipetting
using Serological Pipets!

- Repeat pipetting/multidispensing with less fatigue
- Highly accurate and precise aliquots with one touch
- Now available with Variable aliquoting speed (VS model)
- Use serological pipets but no need to watch the lines
- Plus aspirate and dispense as you're used to

ali-Q is in a class by itself!

Repeat pipetting using serological pipets plus conventional pipetter in one unit



on Variable Speed (VS) models only



set aliquot speed

on all models



set aliquot volume



press aliquot button to dispense each aliquot

