

POWTEQ®



Three-dimensional Mixer M10

◇ Three-dimensional Mixer M10

M10 is suitable for mixing the particulate solid and liquid samples. It can ensure reliable analysis results by rotating evenly, so it is an ideal choice for quality control and R & D.

○ Applications

Chemical materials, semiconductor materials, aviation materials, pharmaceuticals, etc.



M10

○ Features and advantages

- Easy operation , convenient to use
- Low noise and high handling capacity
- Dust-free operation, no environmental pollution

○ Working principle

According to Schatz Geometric theory, the different specific gravity samples are placed into a closed container and kept rolling-over, inverting, rocking in a multi axial space until to the maximum uniformity.

○ Technical data

Speed	20~100rpm	Rated power	200W
Time setting	00 : 01 ~99 : 59 (hr/min)	Power supply	220V , 50/60Hz
Volume	2L/1L/0.5L	Instrument size(W*D*H)	496*640*375mm
Net weight	60kg	Package size(W*D*H)	860*960*780mm