

Zentrifugen von Funke-Gerber
Centrifuges made by Funke-Gerber

FUNKE
GERBER

Nova Safety

Tischzentrifuge zur Fettbestimmung nach Dr. N. Gerber
Bench Centrifuge for the Determination of the Fat Content in Milk
according to Dr. N. Gerber



Nova Safety Gerber Zentrifuge

Eigenschaften:

Deckelverriegelung
 Automatische Bremse (Bremszeit \leftarrow 8 s)
 Timer für Zentrifugierzeit (digital)
 Heizung, auf 65°C geregelt
 Füllmenge: max. 8 Butyrometer

Technische Daten:

RZB: 350 g +/- 50 g
 Drehzahl: 1350 U/min
 Eff. Radius: 170 mm
 Anschluss: 230V/50...60 Hz/500 VA
 Gewicht: 13 kg
 Maße: LxBxH: 470x380x230 mm




Nova Safety Gerber Centrifuge

Features:

Safety lid lock
 Electrical brake time (time \leftarrow 8 s)
 Timer für Zentrifugierzeit (digital)
 Heating, regulated on 65°C
 Suitable for 8 butyrometers (all types)

Technical specifications:

RZB: 350 g +/- 50 g
 Speed: 1350 U/min
 Eff. Radius: 170 mm
 Connection: 230V/50...60 Hz/500 VA
 Weight: 13 kg
 Dimensions: LxBxH: 470x380x230 mm



Bestelldaten/Zubehör:	Art.-Nr.
Zentrifuge Nova Safety	3670

Ordering specifications/Accessories:	Art. No.
Centrifuge Nova Safety	3670

