

Genomic DNA Isolation Reagent Kit

(Blood/Cultured Cell/Tissue)

Catalog Number: SN022-0100

Size: 100 reactions

Lot Number: 40104010

Product Qualification

Specifications:

We show the product pH and recovery ratio in below:

Content	BR	BC	BP
Lot Number	90104011	90104012	90104013
pH	8.1	8.5	5.5
Accredited	√	√	√

Result	Test Sample	Test Volume	Yield (µg)	Agarose check	accredited
1	Whole Blood	200µ l	6.2	√	√
2	Whole Blood	200µ l	6.1	√	√
3	Whole Blood	200µ l	6.5	√	√
4	Whole Blood	200µ l	6.4	√	√
5	Whole Blood	200µ l	6.8	√	√
6	Whole Blood	200µ l	6.8	√	√

Result: Meets specification.

Storage condition: Genomic DNA Isolation Reagent Kit (Blood/Cultured Cell/Tissue) should be stored dry at room temperature (15–25°C).

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses